## TENSIOMÈTRE

FICHE N° 4127

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## Lecomte du Noüy Tensiometer surface and interface

The LECOMTE DU NOÜY Tensiometer is for making all measurements of SURFACE TENSIONS (surface liquid-vapour) and INTERFACE TENSIONS (at the surface of separation of two non-miscible or slightly miscible liquids).

These measurements can be carried out in a few seconds with 2 ml of liquid at about 0.1 dyne.

It can be used for ADSORPTION STUDIES and also lends itself readily to both STATIC and DYNAMIC MEASUREMENTS.

### APPLICATIONS

As follows, we give a summary, necessarily incomplete, of the fields in which this instrument finds applications :

- Scientific Research Laboratories : absolute constants, adsorption, etc.
- Chemical Laboratories : petroleum, oils, rubber, sugar, etc.
- Industrial Laboratories : varnishes, perfumes, soap solutions, volatile liquids, etc.
- Medical and Pharmaceutical Laboratories : physiology, biology, etc.
- Control and Control Laboratories : Navy, Air Force, Army, etc.

### PRINCIPLE OF OPERATION

The LECOMTE DU NOÜY Tensiometer is based on the principle of the torsion balance.

This torsion balance of high sensitivity is intended for measurements of surface and interfacial tensions by the ring method. This method is one of the rare ones which give absolute measurements.

The Tensiometer comprises essentially a torsion thread fixed at one extremity and subjected at the other to a measured torsion; to the thread a bar is fixed and to the bar a ring plunging in the liquid.

The torsion which has to be applied to the thread to restore the ring to its initial reference position after a certain force has been applied to the ring, measures this force.

Readings are obtained directly in dynes per centimeter. Each instrument is adjusted for the length type of ring delivered with it. This adjustment can be refitted readily by the user himself if he changes rings, in conformity with the instructions given in the "mode d'emploi" which accompany each instrument.

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
Figure 1.

## DESCRIPTION

Measurements of surface tension are made by pulling the ring upwards.

On the other hand, interfacial tensions can be done either by upward or downward pulls.

For this purpose, the ring is fixed to the bar *bc* of the parallelogram labelled *abcd* :

the corner *a* of this parallelogram joins the torsion thread ; the corner *d* joins a second thread parallel to the first.

This second thread only acts as a pivot ; only the torsion exerted at *a* produces the pulling couple for the ring. The vertical bar *bc* is detachable and during transportation can be placed in a location provided for this purpose in the transport box for the **Tensiometer**. In operation, it is held by two points (*fig. 2*) fitting into the female cones of the ends of the bars *ab* and *cd*.

The ring is fixed to the hollow rod *e* which extends the bar *bc* under the stirrup arranged to receive the end of the bar *ab* and support it on the cone *b*.

The torsion produced in *a*, in whatever sense, allows a vertical force to be exerted on the ring, either upwards or downwards.

The plane of the ring is always horizontal.

When carrying out absolute measurements, it is necessary that at the time of extraction of the ring, the parallelogram *abcd* should be in the exact position that it was in at the time of calibration.

To obtain this, two fine movements are provided : one which adjusts the height of the surface of the liquid under examination with relation to the **Tensiometer** ; and another which allows the variation of the height of the assembly of the **Tensiometer** with relation to the surface of the liquid.

Surface and interfacial tensions are read on a limb with a double division in dynes per centimeter, on both sides of the zero, and a double vernier at 1/10 (reading directly to 1/10 of a dyne).

Scale and vernier are located in the same plane to avoid any parallax error.

On the vertical rod *bc* is mounted a flag carrying a little lens (*fig. 1*).

A small magnet gives the flag an orientation such that the lens *l* forms the enlarged image of a fixed luminous mark on a ground glass plate forming a central reference mark.

Looking at the dial, the Operator sees this image, and by operating on the control mechanisms of the instrument can maintain this at the height of the reference mark.

This operation can be effected with great precision, since the movements of the arm *bc* cause a displacement of the image of the luminous mark in the same sense and about ten times amplified.

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## THERMOSTATIC CELL

This cell (see cover page) is a double walled receiver of nickel plated brass and is fixed on the plate of the Tensiometer by two screws. With it, water at a known and constant temperature can be circulated around the liquids under examination (the inlet and outlet tubes are not visible in the figure).

It has an aperture for the insertion of a thermometer (not supplied) in the circulation water.

## BIBLIOGRAPHY

A bibliography of the studies carried out with LECOMTE DU NOÛY Tensiometers would comprise a work in itself.

However, the titles of the most important Papers and Memoirs can be furnished on request.

## COMMERCIAL SPECIFICATION

Surface and Interface Tensiometer delivered in working order, without transformer, with a ring, in a fine cabinet : complete ..... SIMPLIN

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### Supplementary Accessories :

Transformer 110-220 V-6 V 50 cycles for the luminous source ..... DAGARBA  
Thermostatic cell ..... THERMO

### Spares :

6 V light source ..... FUNDIUMP  
Platinum ring ..... APLANI  
Average length of this ring : 60 mm (the exact value of the ring average circumference is indicated on the box accompanying each ring).

Certain models of Surface and Interface Tensiometers, previous to that described in the present bulletin, use rings of average length 40 mm. These rings can be supplied on request (Commercial Specification « APLASU »).

They cannot in any circumstances be used with the Tensiometer model « SIMPLIN ».

JOBIN and YVON reserves the right to make in instrument design, any changes which would result suitable in accordance with the latest technical progress.

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